

<Instruction>

Journal : Nature communications

Title : PICASSO allows ultra-multiplexed fluorescence imaging of spatially overlapping proteins without reference spectra measurements

Authors : Junyoung Seo, Yeonbo Sim, Jeewon Kim, Hyunwoo Kim, In Cho, Hoyeon Nam, Young-Gyu Yoon, and Jae-Byum Chang

Affiliation : Korea Advanced Institute of Science and Technology

- **Minimum requirements of the code**

- ✓ To run this code, MATLAB (Mathworks) should be installed. In addition, MATLAB Image Processing Toolbox and Wavelet Toolbox are required to run unmixing code.
- ✓ We have tested the code with MATLAB R2020b running on Windows 10.


- **Demo image**

- ✓ A demo input image is in the 'Data' folder, named '3color_data.tif'.
- ✓ The demo input image was acquired by a single excitation laser but within three different spectral ranges from a mouse brain slice labelled with three preformed rabbit antibody complexes (PV-CF488A ,Neun-ATTO514, GFAP-ATTO532).
- ✓ Simply running the code ('three_color_unmixing.m') will generate an unmixed image, named '3color_data_unmixed.tif' in the 'Results' folder.


- **How to modify the code to unmix other images**

- ✓ Input files should be a tif format images, having three channels where the spectral detection range includes the corresponding fluorophore's emission peak (input *IMG1*, input *IMG 2*, and input *IMG 3* shown in Fig. 3a).

1. **Copy your input image files to the 'Data' folder**



2. **Open the code in the 'Code' folder**



3. **Change the 'filename' as the name of your input image file**

```
imgPath = '../Data/';  
filename = '3color_data.tif';
```

4. Run the unmixing code



5. Resultant unmixed images are saved in the ‘Results’ folder, named ‘[filename]_unmixed.tif’.
6. Input mixed images and unmixed images both can be displayed in individual channel mode and composite mode through the imageJ software program, as shown below.

